

## Plasmid DNA Miniprep Kit

### Description

### Materials Provided:

Material	100 Prep	5 Prep
<b>*Resuspension Buffer</b>	30 ml	1.5 ml (RNase A Solution)
Rnase A Solution (10 µg/µl)	300 µl	1.5 ml
Lysis Buffer	30 ml	1.5 ml
Neutralization Buffer	30 ml	1.5 ml
DNA Binding Buffer (Brown Bottle)	40 ml	2 ml
	20 ml (add 80 ml of ethanol)	1 ml (add 4 ml of ethanol)
<b>**Wash Buffer</b>	96 -100 % ethanol	100 % ethanol
Elution Buffer	20 ml	1 ml
Mini Columns	100	5
Tubes as column inserts	100	5

**Additional requirements:** Microcentrifuge Tubes , 96 - 100% Absolute ethanol.

**Note:** \*\*To Wash Buffer add 96-100% ethanol (In 20 ml of wash buffer + 80 ml ethanol).

**\*Add Rnase A Solution 300 µl to 30 ml Resuspension Buffer.**

### Procedure:

1. Pipette about 1 ml of coli cells into a 1.5 ml microfuge/Eppendorf tubes. Centrifuge the sample at 10,000 rpm for 2-5 minutes at room temperature.
2. Discard the supernatant, and resuspend the cell pellet in 250 µl of Resuspension Buffer containing **RNase A** . Mix by tapping
3. Add 250 µl of Lysis Buffer to the cell **(Do not vortex)**
4. Mix the suspension by gently tapping or by inverting the tube up and down 8-10
5. Add 250 µl of Neutralization Buffer and mix the solution thoroughly by inverting the tube up and down 8-10 **(Do not vortex)**.
6. Centrifuge at 10,000-14,000 rpm for 10 Discard the pellet and save the supernatant.
7. Add 375 µl of DNABinding Buffer to the clear supernatant and
8. Load 550-600 µl of the mixture on to the DNA spin column, centrifuge for 1 -2 minutes and discard the flow

**Note:** You can save the remaining half of the lysate and freeze it at -20°C for

---

**future use. If you plan to use all of it now, this will probably double the amount of the DNA yield.**

9. Wash the DNA spin column with 400  $\mu$ l of Wash Centrifuge the column for 1-2 minutes. Discard the flow through. Wash one more time.

**Date Created**

2024/07/03

LINK BIOTECH