

## Lab consumables

### Description

#### RIA-VIAL and Caps

- RIA VIALS are Poly Propylene tubes generally used in Radioimmuno Assay
- BTL Ria Vials PP are clear and can withstand 3000xg RCF during centrifugation
- The PP tubes can accept common acids alkalies and solvents at room temprature.

#### Products Pack Size Bulk/10case

12x75mm	1x500	500
12x40mm	1x500	350
12x75mm	1x500	600
12x40mm	1x500	450

Previous slide

Next slide

Previous slide

Next slide

#### COLLECTION SWABS

- Sterile dry swabs supplied in shockproof round bottom polypropylene tube.
- Dimensions of tube six inches
- Sterilised by Ethylene oxide
- Expiry date: 48 months from sterilisation date.

#### Products Pack Size Bulk/10case

SWAB TUBE	1x500	3,000
SWAB TUBE	1x500	3,400

#### Sterile Disposable Forceps, PS

- Transparent polystyrene (PS)
- Corrugated handles and rounded serrated tips
- Sterile by Eto Sterilization
- Overall Length (in): 4in

---

**Products   Pack Size**

Forceps, PS 100

[Previous slide](#)[Next slide](#)[Previous slide](#)[Next slide](#)**Sterile Culture Tube with dual position Snap Cap**

Dual position caps provide either aerobic or anaerobic storage of samples. Leave the dual position snap cap in the loose position for aerobic work and in the sealed position for anaerobic cultures.

**Products   Pack Size Bulk/10Case**

Sterile Tube 14ml 500      51272

**Petri Dish**

- High clarity polystyrene
- Flat transparent surface
- Numeric mark on the bottom for easy identification
- Stackable for easy storage and handling

**Products   Pack Size Bulk/10case**

10 Sleeve GR 1 x 480      3,700

20 Sleeve GR 1 x 480      3,600

10 Sleeve E    1 x 480      3,600

20 Sleeve E    1 x 480      3,300

[Previous slide](#)[Next slide](#)[Previous slide](#)[Next slide](#)**CONTACT PLATE (RODAC)**

- Material : Polystyrene Plate : with gridded bottom and protect lid
- Agar bed surface area : 25 cm 2 Gamma irradiated
- Package : 10 per sleeve, 500pcs/case
- Size: 65 X15 mm (outside)

**Products   Pack Size Bulk/10Case**

POLYSTYRENE PLATE 1x500      3,250

## CRYOGENIC STORAGE VIALS

- USP VI Polypropylene
- Autoclavable
- Clear graduations and white marking area
- Non-pyrogenic, DNase/RNase

Products	Pack Size	Bulk/10Case
1.8ml external thread	1x1000	4,500

[Previous slide](#)

[Next slide](#)

[Previous slide](#)

[Next slide](#)

## CRYOGENIC STORAGE BOX

Efficient storage up to temperatures of +90OC, and compatible with a variety of tubes and vials. The cryo cube boxes fit into a variety of Freezer racks, making them versatile. Lock prevents lid opening in storage. Both Microcentrifuge tubes and Cryo vials can be accommodated.

Products	Pack Size	Bulk/10Case
100 PLACES (1.8 ml)	1x4	1,200

## CENTRIFUGE TUBES

Whether centrifugation, mixing, or other applications the BTL Conical Tubes 15 mL and 50 mL will perfectly fit into your workflows.

Products	Pack Size	Bulk/10case
15ml	1x500	3000
50ml	1x500	4800
15ml	1x500	5000
50ml	1x500	5500
50ml	1x500	5500
50ml	1x500	7000

[Previous slide](#)

[Next slide](#)

[Previous slide](#)

Next slide

## Micro Centrifuge Tubes

- Hinged Safe-Lock lid prevents unintentional lid opening during incubation and storage
- Easy labeling on frosted lid and surface on the side
- Functionality is guaranteed from -86 Â°C to 100 Â°C
- Protection of light-sensitive samples is provided by amber colored Safe-Lock Tubes

### Products Pack Size Bulk/10case

0.5ml	1x500	325
1.5ml	1x500	400
2.0ml	1x500	500
1.5ml	1x500	585
2.0ml	1x500	600

## Blood Collection tube

- EDTA K3 provide accurate blood platlet count.
- Efficient anticogulation
- High transparency & visible spray

### Products Pack Size Bulk/10case

EDTA K2	1X100	4800
EDTA K3	1X100	4800
PLAIN CLOT	1X100	4900
FLOURIDE	1X100	4800
S. CITRATE	1X100	6900
GEL	1X100	10500

Previous slide

Next slide

Previous slide

Next slide

## Sputum Container

Sample container made from IP Grade  
99.99% virgin pp use for sample collection, storage  
& Transportation.

### Products Pack Size Bulk/10case

30ml	1X100	175
------	-------	-----

**Products Pack Size Bulk/10case**

30ml sterile	1X100	235
50ml	1X100	290
50ml sterile	1X100	400

**HITACHI CUPS**

The nesting sample cup is used when smaller sampling is necessary. It is used in conjunction with other test tubes or original blood collection tubes. To use, simply transfer the sample from the original collection tube into the nesting cup. Then, place the nesting cup inside the original collection tube. The nesting cup “rides” along with the originally labeled/barcoded tube in the analyzer. This procedure saves time by eliminating the need to re-label the smaller sample.

**Products Pack Size Bulk/10case**

2ml	1x500	500
-----	-------	-----

Previous slide

Next slide

Previous slide

Next slide

**MULTIPURPOSE VIALS**

Multi Purpose Vial. Made up of Polypropylene. Specially designed for transportation of sample / serum.

**Products Pack Size Bulk/10case**

5ml	1x500	700
-----	-------	-----

**Cryovials**

- Pack Size: 500 Pcs
- Material: PP
- Internally Threaded Cap Style
- Large White Frosted portion for easy writing

Previous slide

Next slide

Previous slide

Next slide

## Pipettes

- Ergonomic Design
- Digital Display
- Calibration can be done easily with tool supplied
- Light Weight
- Autoclavable

## Micro Cover Glass

Small Squares of glass that cover the specimen are placed on the Microscope Glass Slides & are packed in controlled conditions. Packaging- Approx 10gm of transparent coverslips are packed in 1 packet and such 20 packets are placed in 1 box.

Previous slide

Next slide

Previous slide

Next slide

## DOUBLE MICROSCOPE GLASS SLIDES

- Packaging- Double Packet (Top & Bottom Separate)
- Size:- 76.2mm x 25.4mm (1"x3")
- Thickness:- 1mm-1.2mm
- 1 packet contains 50 pieces(Nos.) of slides and such 50 packets are packed in a corrugated box.

## FROSTED MICROSCOPE SLIDES

- Frosted at One End – Both Sides
- Size:- 76.2mm x 25.4mm (1"x3")
- Thickness:- 1mm-1.2mm
- 1 packet contains 50 pieces(Nos.) of slides and such 50 packets are packed in a corrugated box.

Previous slide

Next slide

Previous slide

Next slide

## PLAZA DELUXE MICROSCOPIC SLIDES

- Ground edges and lint-free packing
- Size:- 76mm x 26mm (+/-1mm)
- Thickness:- 1.25mm (+/-0.1mm)
- 1 packet contains 50 pieces(Nos.) of slides and such 50 packets are packed in a corrugated box.

## Date Created

2024/07/24