

anti- TAK1 antibody

Description

Product Name: anti- TAK1 antibody

Catalogue No.: FNab08492

Size: 100 µg

Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Clonality: polyclonal

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

IMMUNOGEN INFORMATION


Immunogen: mitogen-activated protein kinase kinase kinase 7

Observed MW: 64-67 kDa, 75-80 kDa


Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

IMAGES



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using BS-29480BANF (MAP3K7 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



CHO cells were subjected to SDS PAGE followed by western blot with BS-29480BANF(MAP3K7 antibody) at dilution of 1:1000

LINK BIOTECH

Date Created
2024/07/02