


15 Lipoxygenase 1 Antibody

Description

Product Name	15 Lipoxygenase 1 Antibody
Catalogue No.	DF13494
Size	100ul
Form	liquid
Purification	The antiserum was purified by peptide affinity chromatography using Sulfolink Coupling Resin (Thermo Fisher Scientific).
Source:	Rabbit
Clonality	Polyclonal
Specificity:	15 Lipoxygenase 1 Antibody detects endogenous levels of total 15 Lipoxygenase 1.
Conjugate:	Unconjugated.
Storage	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 Â°C. Stable for 12 months from date of receipt.

Western blot analysis of extracts from Rat liver, using 15 Lipoxygenase 1 Ab. The lane on the left was treated with blocking peptide.



Date Created

2024/07/01