

Goat Anti-Rabbit IgG (H+L) Fluor488-conjugated

Description

Product Name	Goat Anti-Rabbit IgG (H+L) Fluor488-conjugated
Catalogue No.	S0017
Size	100ul
Form	liquid
Purification	NA
Source:	Goat
Clonality	Polyclonal
Specificity:	
Conjugate:	Fluor488. Ex:499nm/Em:520nm
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. Store at -20â„ƒ. Avoid repeated thawing and freezing.

Date Created

2024/07/03

LINK BIOTECH