

ELISA Kit for Interleukin 6 (IL6)

Description

| | |
|----------------------|---|
| Catalog No. | SEA079Mu |
| Unit Size | 48T, 96T, 96T*5, 96T*10, 96T*100 |
| Product Category | ELISA (Quantitative) |
| Reactivity | Mus musculus (Mouse) |
| Sample | serum, plasma, tissue homogenates, cell lysates, cell culture supernates biological fluids |
| Assay Method | Double-antibody Sandwich |
| Assay Duration | 3h |
| Sensitivity | The minimum detectable dose of this kit is typically less than 2.7pg/mL. |
| Standard Curve Range | 7.8-500pg/mL |
| Storage/Storage | -20°C/1 year; 4°C/6 months |

Date Created

2024/07/03

LINK BIOTECH