

BSiFuge – D06

Description

SPECIFICATIONS

BSIFUGE D06

Fixed Angle Rotor (Optional Rotors)	6Å—15 ml(Autoclavable) – 3684 g 8Å—15 ml(Autoclavable – 3873 g
Max RPM	6500
Motor Type	Brushless DC
Speed Setting	500 to 6500 RPM
Timer Setting	1 to 30 mins & infinite mode
Lid Lock Safety	Yes
Noise Level	< 60 dB
Imbalance Detection	Yes
Dimensions (W x D x H)	260 x 244 x 205 mm
Weight	4 kg
Adaptors (included in standard pack)	Chosen rotor with (adaptors for 5-7 ml & 3/4 ml tubes – takes most standard blood collection tubes)

Date Created

2024/07/04