

BSiRupt 24P

Description

Agitation speed range	0.80 m/s ~ 7.50 m/s (in steps of 0.05 m/s) 490 rpm â€“ 4550 rpm (in step of 30 rpm)
Run Time	1 second to 90 second (in 1 second increment)
Pause Time	10 second to 600 second (in 1 second increment)
Number of Samples (Capacity)	1~24 samples (24 x 2 ml, 12 x 7 ml)
Number of Cycles	10
Programmable (Stored)	99 Saved Programs
Dimensions (W x D x H)	350 x 490 x 315 mm
Weight	32.5 kg
Power Supply	230 VAC/ 50Hz
Power Consumption	160 W

Date Created

2024/07/04

LINK BIOTECH