

BSiShak BL Uno VT

Description

Speed Setting	300 to 1800 RPM
Orbital Diameter	2 mm
Motor type	BLDC
Timer setting	1 to 180 mins and infinity
Digital Display	Yes
Inbuilt Counterbalance	Yes
Dimension (W x D x H)	127 x 86 x 85 mm
Weight	1.36 kg
Attachments (included in standard park)	Replaceable clamps to use different types of standard micro plates
Included Attachment	Attachment for 20 x 1.5 / 2.0 ml micro tubes

Date Created

2024/07/04

LINK BIOTECH