

## BSiRoll DR24

### Description

Capacity (supplied with any one chosen disc plate)	12 tubes x 50 ml 24 tubes x 15 ml 44 tubes x 2 ml
Programmable Mode	Yes. User can set and save upto 9 User-Defined programs
Variable Speed	10-80 RPM
Run Time	1-99 mins & infinite
Motor Type	BLDC motor
Pulse Mode	Yes. Alternating clockwise and anticlockwise motion with 2 sec stops for better and thorough mixing.
Angle of Disc Plate	Adjustable from 0-90°
Max Load (Capacity)	Upto 2 kg
Digital Display	Yes
Input Voltage	24V
Dimension	383 x 325 x 325 mm
Protection class	IP 21

### Date Created

2024/07/04