

BSiShak 3D-5 NXT

Description

Weight Capacity	5 to 7 kg
Variable Speed	15-120 RPM
Movement Type	Tumbling
Tumbling Angle	7 ⁰
Timer Range	1-999 mins & infinite mode
Operation Type	Timed & Continuous Operation
Pulse Mode	Yes (Run-Pause-Run Mode)
Plate Size	350 mm x 350 mm
Power Consumption	24W
Attachments	Spiked Rubber Mat (Included) Sticky Pad Mat (Optional)

Date Created

2024/07/04

LINK BIOTECH