

Pipette Controller



Description

TECHNICAL INFORMATION

Power	2 x 1.2 V / 1000 mAh rechargeable NiMH batteries.
Battery usage Duration	Eight hours continuous use
Battery Charging Duration	2 – 3 hours to charge fully
Volumetric Range	Plastic or Glass pipette from 0.1 ml up to 200 ml
Safety System	Replaceable hydrophobic 0.45 µm PTFE filter & Safety check valve


Standard Accessories	<ul style="list-style-type: none"> NiMH Batteries Battery Charger Hydrophobic PTFE membrane Filter 0.45 µm – 2Pc Hydrophobic PTFE membrane Filter 0.20 µm – 1Pc Wall mounting stand Desktop stand
----------------------	---

ORDERING INFORMATION

DESCRIPTION	CAT#	COLO
Pipet Controller – Red	PH01-R	
Pipet Controller – Blue	PH01-B	


Pipet Controller - Orange

PH01-O




Pipet Controller - Green


PH01-G

PTFE Membrane Filter Set, 0.20 μm (5 pcs)

PH01-12


PTFE Membrane Filter Set, 0.45 μm (5 pcs)

PH01-13




Battery Charger

PH01-14




Silicone Pipette adaptor with safety valve

PH01-10




Desktop stand

PH01-18




Wall mounting stand

PH01-17



NiMH battery (set of two batteries)

PH01-16



Date Created
2024/07/03

LINK BIOTECH