

## BSiSwix

### Description

Max Speed	2800 RPM
Orbital diameter	4 mm
Motor type	AC Motor
Max Load Capacity	500 gm
Speed accuracy	Â± 5%
Different modes Touch, Stand-by & Continuous mode	Touch, Stand-by & Continuous mode
Dimensions	205 x 136 x 138.5 mm
Power Supply	220 â€“ 240VAC, 50Hz
Ambient temperature	5 to 40 <sup>0</sup> C
Weight	3 kg
Optional Attachments (not included)	Universal platform for attaching either microplate or microtube or flask Attachments

### Date Created

2024/07/04

LINK BIOTECH