

0.5 ml Micro Centrifuge Tubes

Description

General Features

- Manufactured from FDA approved resin
- Certified RNase, DNase, Human DNA, Endotoxin & Pyrogen free
- Can be centrifuged at up to 20,000 RCF
- Withstand temperatures from -80 to +121C
- Fit all standard and high-capacity centrifuge rotors
- Supplied in double safety self-standing bulk bags

Cat. No.	Description
BCT-005-B	0.5 ml Micro Centrifuge Tube, Natural
BCT-005-B-S	0.5 ml Micro Centrifuge Tube, Natural, Sterile
BCT-005-L-B-S	0.5 ml Micro Centrifuge Tube, Low Retention, Sterile
BCT-005-B-C	0.5 ml Micro Centrifuge Tube, Color (C= GREEN/ ORANGE/ PINK/ RED/ YELLOW/ VIOLET)
BCT-005-B-S-C	0.5 ml Micro Centrifuge Tube, Color, Sterile (C= GREEN/ ORANGE/ PINK/ RED/ YELLOW/ VIOLET)
BCT-005-B-A	0.5 ml Micro Centrifuge Tube,Amber color
BCT-005-B-S-A	0.5 ml Micro Centrifuge Tube,Amber color, Sterile

Date Created

2024/07/03