

Actidione, for tissue culture, 98%+

Description

DESCRIPTION

Cycloheximide (Actidione) is a eukaryote protein synthesis inhibitor, produced by the bacterium Streptomyces griseus.

APPLICATION

Actidione, for tissue culture, 98%+ Cas 66-81-9 is used only in in vitro research applications, and is not suitable for human use as a therapeutic compound.

PROPERTIES

Categories

Assay  
Appearance (Colour)  
Appearance (Form)  
Solubility  
Melting point

Antibiotics

98%+  
White to off white  
Powder  
Ethanol: soluble, clear to hazy  
110-119 °C

Code

6141A-100  
6141A-1G  
6141A-250G  
6141A-C

Packings

100 mg  
1 gm  
250 gm  
Custom size

Date Created

2024/07/04