

Actidione, for tissue culture, 98%+

Description

DESCRIPTION

Cycloheximide (Actidione) is a eukaryote protein synthesis inhibitor, produced by the bacterium Streptomyces griseus.

APPLICATION

Actidione, for tissue culture, 98%+ Cas 66-81-9 is used only in in vitro research applications, and is not suitable for human use as a therapeutic compound.

PROPERTIES

Categories

Assay
Appearance (Colour)
Appearance (Form)
Solubility
Melting point

Antibiotics

98%+
White to off white
Powder
Ethanol: soluble, clear to hazy
110-119 °C

Code

6141A-100
6141A-1G
6141A-250G
6141A-C

Packings

100 mg
1 gm
250 gm
Custom size

Date Created

2024/07/04