



ANTI-RABBIT SERUM (GOAT)

Description

PRODUCT NAME: ANTI-RABBIT SERUM (GOAT)

PRODUCT NUMBER: 200840

BATCH NUMBER: 65857

FORM: Liquid (Polyclonal, Delipidified and 0.2micron filtered)

BUFFER COMPOSITION: 20mM Tris at pH 7.4 containing 250mM NaCl and 15mM Na₃

MOLECULAR WEIGHT: 84.31

STORAGE: 2°C-8°C

MFG. DATE: MAR.2023

EXP. DATE: APR.2027

SOURCE ANIMAL HEALTH: Inspected and Certified free of Rinderpest, Foot mouth disease,

Contagious bovine pleuropneumonia and negative for blue tongue.

Test	Test Method	Acceptance Criteria
Titre	Ouchterlony	≥ 1:8
Total Protein	Biuret Method	50-150 mg/ml
Specificity	Immuno Electrophoresis	Immuno Electrophoresis gives precipitin are specific when re Normal rabbit serum.

Caution : To be used for in-vitro diagnostic purpose only:

Date Created
2024/07/02