

## ANTI-RABBIT SERUM (GOAT)

### Description

**PRODUCT NAME:** ANTI-RABBIT SERUM (GOAT)

**PRODUCT NUMBER:** 200840

**BATCH NUMBER:** 65857

**FORM:** Liquid(Polyclonal, Delipidified and 0.2micron filtered)

**BUFFER COMPOSITION:**

20mM Tris at pH 7.4 containing 250mM NaCl and 15mM NaN3

**MOLECULAR WEIGHT:** 84.31

**STORAGE:** 2-8C

**MFG. DATE:** MAR.2023

**EXP.**

**DATE:** APR.2027

**SOURCE ANIMAL HEALTH:** Inspected and Certified free of Rinderpest, Foot mouth disease,

Contagious bovine pleuropneumonia and negative for blue tongue.

Test	Test Method	Acceptance Criteria
Titre	Ouchterlony	> 1:8
Total Protein	Biuret Method	50-150 mg/ml
Specificity	Immuno Electrophoresis	Immuno Electrophoresis gives precipitin arcs specific when compared to normal rabbit serum.

**Caution:**

To be used for in-vitro diagnostic purpose only:

**Date Created**

2024/07/02