

## anti- HBEGF antibody

### Description

|               |   |
|---------------|---|
| Product Name  | anti- HBEGF antibody  |
| Catalogue No. | 12601AB   |
| Size          | 100 µg  |
| Form          | liquid  |
| Purification  | Immunogen affinity purified   |
| Purity        | ≥95% as determined by SDS-PAGE  |
| Host          | Rabbit  |
| Clonality     | Polyclonal  |
| Isotype       | IgG   |
| Storage       | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### IMMUNOGEN INFORMATION


|             |   |
|-------------|---|
| Immunogen   | Proheparin-binding EGF-like growth factor |
| Observed MW | 23 kDa                                    |

### APPLICATION


|                      |                                   |
|----------------------|-----------------------------------|
| Tested Application   | ELISA, WB, IHC                    |
| Recommended Dilution | WB: 1:500-1:2000; IHC: 1:50-1:200 |

### IMAGES

Immunohistochemistry of paraffin-embedded Human colon cancer using (HBEGF antibody) at dilution of 1:1000



293T cell lysates were subjected to SDS PAGE followed by western blot with (HBEGF antibody) at dilution of 1:1000



**Date Created**  
2024/07/01

LINK BIOTECH