

anti- Surfactant protein D antibody

Description

Product Name	anti- Surfactant protein D antibody
Catalogue No.	10480BSAB
Size	100µg
Form	liquid
Purification	Immunogen affinity purified
Purity	≥95% as determined by SDS-PAGE
Host	Rabbit
Clonality	polyclonal
Isotype	IgG
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

IMMUNOGEN INFORMATION


Immunogen	surfactant protein D
Observed MW	37 kDa

APPLICATION


Tested Application	ELISA, IHC, IP
Observed MW	37 kDa

IMAGES

IP Result of anti-Surfactantprotein D (IP: 10480BSAB,3ug; Detection: 10480BSAB1:300) with mouse lung tissuelysate 2800ug.



Immunohistochemical of paraffin-embedded human lung cancer using 10480BSAB (Surfactant protein D antibody) at dilution of 1:50



Date Created
2024/07/02