

## anti- TAK1 antibody

### Description

**Product Name:** anti- TAK1 antibody

**Catalogue No.:** FNab08492

**Size:** 100 µg

**Form:** liquid

**Purification:** Immunogen affinity purified

**Purity:** ≥95% as determined by SDS-PAGE

**Clonality:** polyclonal

**Storage:** PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### IMMUNOGEN INFORMATION


**Immunogen:** mitogen-activated protein kinase kinase kinase 7

**Observed MW:** 64-67 kDa, 75-80 kDa


**Tested Application:** ELISA, IHC, WB

**Recommended dilution:** WB: 1:500-1:2000; IHC: 1:20-1:200

### IMAGES



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using BS-29480BANF (MAP3K7 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



CHO cells were subjected to SDS PAGE followed by western blot with BS-29480BANF(MAP3K7 antibody) at dilution of 1:1000

LINK BIOTECH

**Date Created**  
2024/07/02