

## Biosimilar PK ELISA Kits

### Description

CATALOGUE NO	PRODUCT	ASSAY MATRIX	ASSAY RANGE	SENSITIVITY	STORAGE
PK1001	<b>Infliximab (Remicade®) PK ELISA KIT</b>	Human serum and plasma	6.58-50 ng/ml	6 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1002	<b>Adalimumab (Humira®) PK ELISA KIT</b>	Human serum and plasma	5.2 - 200 ng/ml	5.2 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1003	<b>Etanercept (Enbrel®) PK ELISA KIT</b>	Human serum and plasma	300 ng/ml to 2.34 ng/ml	2.34 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1004	<b>Rituximab (Mabthera®) PK ELISA KIT</b>	Human serum and plasma	6.25-200 ng/ml	6.25 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1005	<b>Bevacizumab (Avastin®) PK ELISA KIT</b>	Human serum and plasma	0.1 - 3.1 ug/ml	1.17 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1006	<b>Ranibizumab (Lucentis®) PK ELISA KIT</b>	Human serum and plasma	1250 - 9.77 ng/ml	4.88ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C
PK1009	<b>Denosumab (Xgeva®) PK ELISA KIT</b>	Human serum and plasma	1.95 - 1000 ng/ml	3.91 ng/ml	Shipped at 2-8°C. Longterm Storage at -20°C

CATALOGUE NO	PRODUCT	ASSAY MATRIX	ASSAY RANGE	SENSITIVITY	STORAGE
IM5004	Rituximab (Mabthera®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cutoff values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C
IM5005	Adalimumab (Humira®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cut off values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C
IM5006	Ranibizumab (Lucentis®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cut off values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C

Date Created  
2024/07/03

LINK BIOTECH