

Biosimilar PK ELISA Kits

Description

CATALOGUE NO	PRODUCT	ASSAY MATRIX	ASSAY RANGE	SENSITIVITY	STORAGE
PK1001	Infliximab (Remicade®) PK ELISA KIT	Human serum and plasma	6.58-50 ng/ml	6 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1002	Adalimumab (Humira®) PK ELISA KIT	Human serum and plasma	5.2 â€“ 200 ng/ml	5.2 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1003	Etanercept (Enbrel®) PK ELISA KIT	Human serum and plasma	300 ng/ml to 2.34 ng/ml	2.34 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1004	Rituximab (Mabthera®) PK ELISA KIT	Human serum and plasma	6.25-200 ng/ml	6.25 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1005	Bevacizumab (Avastin®) PK ELISA KIT	Human serum and plasma	Â 0.1 â€“ 3.1 ug/ml	1.17 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1006	Ranibizumab (Lucentis®) PK ELISA KIT	Human serum and plasma	1250 â€“ 9.77 ng/ml	4.88ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C
PK1009	Denosumab (Xgeva®) PK ELISA KIT	Human serum and plasma	1.95 â€“ 1000 ng/ml	3.91 ng/ml	Shipped at 2-8Â°C. Longterm Storage at -20Â° C

CATALOGUE NO	PRODUCT	ASSAY MATRIX	ASSAY RANGE	SENSITIVITY	STORAGE
IM5004	Rituximab (Mabthera®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cutoff values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C
IM5005	Adalimumab (Humira®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cut off values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C
IM5006	Ranibizumab (Lucentis®) Immunogenicity PK ELISA KIT	Human serum and plasma	refer cut off values	refer cut off values	Shipped at 2-8°C. Longterm Storage at -20°C

Date Created
2024/07/03

LINK BIOTECH