

## BSiFuge C4000 NXT

### Description

Dimension (W x D x H)	381 x 480 x 185 (mm)
Weight of Equipment	8.2 kg (without rotor)
Permissible Temp.	5-40°C
Humidity	<80%
Altitude	Use upto an altitude of 2000 m above MSL
Pollution Degree	Â 2
Environment	For indoor use only
Noise Level	â‰‰65 dB
Power Adaptor Details	Input: 100-240 VAC, 50/60 Hz Output: 24V 5A
Power Consumption	90 W

	<b>ROTOR SPECIFICATION</b>	<b>6 TUBES ROTOR</b>	<b>8 TUBES ROTOR</b>	<b>16 TUBES ROTOR</b>
Rotor		6 x 10ml Swing Rotor	8 x 15ml Fix angle Rotor	16 x 10ml Fix angle Rotor
Maximum Capacity		6 x 2ml, 6 x 3ml, 6 x 4ml 6 x 10ml, 6 x 9ml, 6 x 6ml	8 x 15ml Conical tube	16 x 10ml Gel tube
Nominal RCF		1840 gÂ   2270 g	1950 g	2180 g
Acceleration Time		45 second	45 second	45 second
Deceleration Time		40 second	70 second	70 second
Power		90 W	24 W	35 W
Current Consumption		3.7 A	1.01 A	1.45 A

### Date Created

2024/07/04