

## BSiFuge CS03 NXT

### Description

Capacity	2 x 50 ml, 3 x 50 ml (swing out rotor)
Max RPM/RCF (Xg)	3200 RPM / 1470g
Motor Type	Brushless DC motor
Digital Display	Yes
Speed Setting	500-3200 RPM
Timer Setting	1-999 mins & infinite mode
Programmable Feature	User can set and save upto 99 user-defined programs (protocols)
Lid Lock Safety & Imbalance Detection	Yes
Dimensions (L x W x H)	381 x 480 x 185 mm
Weight	8.2 kg (without rotor)
Voltage / Consumption	24 V 1.5A / 36 W
<b>Date Created</b>	
2024/07/04	

LINK BIOTECH