

BSiFuge CS2P NXT

Description

Capacity	2 micro plate (96 & 384 wells)
Max RPM/RCF (Xg)	2200 RPM/560g
Motor Type	Brushless DC motor
Digital Display	Yes
Speed Setting	500-2200 RPM
Timer Setting	1-999 mins & infinite mode
Programmable Feature	User can set and save upto 99 user-defined programs (protocols)
Lid Lock Safety & Imbalance Detection	Yes
Dimension (L x W x H)	381 x 480 x 185 mm
Voltage / Consumption	24V – 1.4A / 35W
Weight	8.2 kg (without rotor)
Date Created	
2024/07/04	

LINK BIOTECH