

BSiFuge CYTO

Description

Dimension (W x D x H)	381 x 480 x 185 mm
Weight of Equipment	8.3 kg (without rotor)
Permissible Temp.	5-40°C
Humidity	<80%
Altitude	Use upto an altitude of 2200 m above MSL
Pollution Degree	2
Environment	For indoor use only
Noise Level	≤65 dB
Power Adaptor Details	Input: 100-240 VAC, 50/60 Hz Output 24V 5A
Power Consumption	48 W

Date Created

2024/07/04

LINK BIOTECH