

## BSiRoll PR35

### Description

Capacity (supplied with any one chosen disc plate)	16 tubes x 50 ml 35 tubes x 15 ml 55 tubes x 2 ml
Programmable Mode	Yes. User can set and save up to 9 User-Defined programs
Variable Speed	10-80 RPM
Run Time	1-99 mins & infinite
Motor Type	BLDC motor
Pulse Mode	Yes. Alternating clockwise and anticlockwise motion with 2 sec stops for better and thorough mixing.
Angle of Disc Plate	Adjustable from 0-90 degrees
Weight (Capacity)	1.5 kg
Digital Display	Yes
Input Voltage	24V
Dimension (W x D x H)	540 x 274 x 186 mm
Protection class	IP 21

### Date Created

2024/07/04