

BSiRoll RnR4

Description

Capacity	9 x 3/4 ml tubes
Motor Type	AC Synchronous
Speed	30 RPM
Run Time	Continuous / Infinite
Tilt angle	+/- 3°
Speed Accuracy	+/- 2 RPM
Noise Level	Noiseless Operation
Load Capacity	1.5 kg
Roller Size (mm)	300 (length) x 30 (roller dia) mm
Protection class	IP 21
Dimension (W x D x H)	475 x 210 x 80 mm
Weight	2 kg

Date Created

2024/07/04

LINK BIOTECH