

BSiShak PS 10/20

Description

| | | | | |
|-------------------------|--|---------|---------------|---------|
| Shaking Movement | Linear | Orbital | Linear | Orbital |
| Orbital Diameter | 10 mm | | 20 mm | |
| Max Weight | 7.5 kg | | | |
| Speed | 50 to 500 RPM | | 50 to 300 RPM | |
| Motor Type | Brushless DC Motor | | | |
| Time Setting Range | Upto 99:59 hrs | | | |
| Platform Size | 330 x 330 mm | | | |
| Programmable Feature | 6 programs with each program having 9 sub-programs | | | |
| Program Looping | Yes | | | |
| Dimensions of Equipment | 408 x 350 x 99 mm | | | |

Date Created

2024/07/04

LINK BIOTECH