

BSiShak Quattro Nxt

Description

Motor Type	BLDC motor
Type Of Motor	Orbital
Orbital Diameter	3 mm
Plate Capacity	4 Micro Plate
Speed Setting	200 to 1200 RPM
Timer setting	1 to 999 mins & infinite mode
Pulse Mode	Yes (clockwise & counter clockwise)
Dimension (W x D x H) mm	267x361x155
Weight	5.4 kg
Inbuilt Counterbalance	Yes
Display	LED

Date Created

2024/07/04

LINK BIOTECH