

BSiShak TS4 NXT

Description

Motor Type	BLDC motor
Type of Movement	Orbital
Orbital Diameter	3 mm
Number of plates	Up to 4
Temperature Range	Ambient +3 to 70°C
Temperature Stability	±0.5°C
Temperature Accuracy	±0.5°C
Speed Setting	200 to 1200 RPM at every 10 RPM
Time Setting	1 to 99 h 59 min & infinity
Pulse mode	Yes (clockwise & counter clockwise)
Display	LED
Dimension (W x D x H) mm	267 x 361 x 155
Permissibly Ambient Temperature	4 to 40°C
Weight	8.45 kg

Date Created

2024/07/04