

BSiShak US15

Description

Shaking Movement	Orbital / Linear
Orbital Diameter	30 , 40 , 50, 60 & 70 mm
Max Operational Weight	15 kg
Speed	50-300 RPM (30/40/50 mm) & 50-250 RPM (60/70 mm)
Motor Type	Brushless DC Motor
Time Setting Range	Up to 99 hours & 59 minutes
Platform Size	745 x 490 mm
Programmable Feature	6 programs with each program having 9 sub-programs
Program Looping	Yes
Dimensions (W x D x H)	825 x 625 x 165 mm

Date Created

2024/07/04

LINK BIOTECH