

BSiStir 200L

Description

Motor Type	Brushless DC (BLDC)
Max. Capacity	200 litres
Variable Speed	100 â€“ 600 RPM
Run Time	1-99:59 mins & infinite
Ramp up Speed	User can select different Acceleration rates (A1-A9 with A9 being the fastest rate and A1 slowest)
Pulse Mode	Yes. Alternating clockwise and anticlockwise motion with 3 sec stop for better and thorough mixing.
Digital Display	Yes
Accuracy	Â± 10 RPM
Plate Dimensions (SS)	250 x 250 mm
Input Voltage	240 V
Max. Stir bar length	120 mm
Protection class	IP 21
Structure frame sizeÂ (With adjustable height)	505 x 500 x 105 mm

Date Created

2024/07/04

LINK BIOTECH