

BSiStir 50L

Description

Motor Type	Brushless
Max. Capacity	50 litres
Variable Speed	100 – 1000 RPM
Run Time	1-99:59 mins & infinite
Ramp up Speed	User can select different Acceleration rates (A1-A9 with A9 being the fastest rate and A1 slowest)
Pulse Mode	Yes. Alternating clockwise and anticlockwise motion with 3 sec stop for better and thorough mixing.
Digital Display	Yes
Accuracy	± 10 RPM
Plate Dimensions (SS)	184 x 184 mm
Input Voltage	240 V
Max. Stir bar length	80 mm
Protection class	IP 21
Structure frame size (With adjustable height)	394 x 392 x 101 mm

Date Created

2024/07/04