

BSiStir HP 10M

Description

Motor Type	BLDC
Number of Stirring Positions	10
Stirring Quantity Max	400 ml per position
Variable Speed	300-1500 RPM in steps of 10 RPM
Heater Temperature Range	Ambient to 120Â° C
Run Time	1-999 mins & infinite
Pulse Mode	Programmable pulse mode â€” user can set time from 30 to 99 secs
Set-up Plate Material	Stainless steel plate with silicone sheet
Set-up Plate Dimensions	460 x 178 mm
Weight	5 kg
Protection Class	IP 42

Date Created

2024/07/04

LINK BIOTECH