

## BSiStir HP 5M

### Description

|                              |                                                                |
|------------------------------|----------------------------------------------------------------|
| Motor Type                   | BLDC                                                           |
| Number of Stirring Positions | 5                                                              |
| Stirring Quantity Max        | 400 ml per position                                            |
| Variable Speed               | 300-1500 RPM in steps of 10 RPM                                |
| Heater Temperature Range     | Ambient to 120° C                                              |
| Run Time                     | 1-999 mins & infinite                                          |
| Pulse Mode                   | Programmable pulse mode - user can set time from 30 to 99 secs |
| Set-up Plate Material        | Stainless steel plate with silicone sheet                      |
| Set-up Plate Dimensions      | 456 x 118 mm                                                   |
| Weight                       | 3 kg                                                           |
| Protection Class             | IP 42                                                          |

### Date Created

2024/07/04

LINK BIOTECH