

BSiSwix Jr. NXT and BSiSwix Jr. VT NXT

Description

|                                 |                        |
|---------------------------------|------------------------|
| Motor Type                      | DC                     |
| Type Of Movement                | Orbital                |
| Maximum Speed                   | 4200RPM   2000-4200RPM |
| Shaker Diameter                 | 4.5 mm                 |
| Tube Diameter                   | Up to 30 mm            |
| Speed Accuracy                  | ± 250 RPM              |
| Operation Type                  | Touch                  |
| Protection Class                | IP 43                  |
| Permissible Ambient Temperature | 4°C to 45 °C           |
| Weight                          | 0.75 Kg                |
| Permissible Relative Moisture   | <80%                   |
| Noise Level                     | <50 dB                 |
| Dimensions (L X W X H)          | 120 x 120 x 103mm      |
| Input Voltage                   | 100-240VAC 50/60Hz     |
| Power Consumption               | 52W                    |

Date Created  
2024/07/04