

## BSiSwix MV

### Description

Attachments	<ul style="list-style-type: none"><li>• 7 ml tube foam to accommodate 60 tubes</li><li>• 15 ml tube foam to accommodate 45 tubes</li><li>• 50 ml tube foam to accommodate 18 tubes</li><li>• 250 ml flask foam to accommodate 6 flasks</li><li>• 1.5/2 ml tube foam to accommodate 72 tubes</li><li>• 50 ml tube foam (quechers assay) to accommodate 7 tubes</li></ul>
Digital Display	Yes
Variable Speed	300-2500 RPM
Motor Type	BLDC
Run Time	1-999 mins & infinite
Motion Type	Clockwise (CW) and Counterclockwise (CCW) motion
Pulse Mode	Programmable Pulse with selection of variable time between 30-99 seconds
Weight	18 kg
Input Voltage	24 V
Power Consumption	75 W
Protection class	IP 21
<b>Date Created</b>	
2024/07/04	