

BSiSwix VT

Description

|                                     |                                                                                       |
|-------------------------------------|---------------------------------------------------------------------------------------|
| Max Speed                           | 4200 RPM                                                                              |
| Orbital diameter                    | 4 mm                                                                                  |
| Motor Type                          | Brushless DC                                                                          |
| Digital Display                     | Yes                                                                                   |
| Dimensions (W X D X H)              | 205 x 136 x 138.5 mm                                                                  |
| Speed Setting                       | 300-4200 RPM                                                                          |
| Timer Setting                       | 1-999 mins (min/sec display)                                                          |
| Pulse mode                          | Yes, Programmable Pulse                                                               |
| Max Load Capacity                   | 500 gm                                                                                |
| Weight                              | 2.8 kg                                                                                |
| Optional Attachments (not included) | Universal platform for attaching either 96 microplate / microtube / flask Attachments |

Date Created

2024/07/04

LINK BIOTECH