

Cell Culture Dish

Description

Cell culture dishes are made from polystyrene (PS) materials and are suitable for laboratory inoculation, marking, weighing, separation and treatment of tissue or cells.

PRODUCT FEATURE

- â– Uniform thickness, no distortion for excellent observation.
- â– Easy to grip with gear ring design, reduce contamination risk.
- â– The ring protrusion on the cover is closely combined with the bottom to facilitate storage and reduce medium volatilization.
- â– Open petri dish lids with regular protrusions are available.
- â– Vacuum plasma surface treatment for ideal adhesive performance.
- â– No pyrogen, endotoxin free, no cytotoxic.
- â– Sterile by Gamma irradiation, SAL<10⁻⁶.

PRODUCT PARAMETER

Item No.	Size(mm)	Min package	Sale package
C070035	35	10pcs/pack	50packs/case
C070060	60	10pcs/pack	50packs/case
C070100	100	10pcs/pack	30packs/case
C070150	150	5pcs/pack	20packs/case

Date Created

2024/07/03