

## Cell Culture Dish

### Description

Cell culture dishes are made from polystyrene (PS) materials and are suitable for laboratory inoculation, marking, weighing, separation and treatment of tissue or cells.

### PRODUCT FEATURE

- Uniform thickness, no distortion for excellent observation.
- Easy to grip with gear ring design, reduce contamination risk.
- The ring protrusion on the cover is closely combined with the bottom to facilitate storage and reduce medium volatilization.
- Open petri dish lids with regular protrusions are available.
- Vacuum plasma surface treatment for ideal adhesive performance.
- No pyrogen, endotoxin free, no cytotoxic.
- Sterile by Gamma irradiation, SAL $\leq 10^{-6}$ .

### PRODUCT PARAMETER

Item No.	Size(mm)	Min package	Sale package
C070035	35	10pcs/pack	50packs/case
C070060	60	10pcs/pack	50packs/case
C070100	100	10pcs/pack	30packs/case
C070150	150	5pcs/pack	20packs/case

### Date Created

2024/07/03