

## Copper Microplate Assay Kit BSK1154

### Description

#### PRODUCT INFORMATION:-

☛ **Product Name** :-Copper Microplate Assay Kit

☛ **Catalogue No.** :- BSK1154

☛ **Source** :-

☛ **Reactivity** :- All

☛ **Application** :- Functional Assay

☛ **Description** :- Detection and Quantification of Glutathione Content.

☛ **Alternative Names** :- Glutathione, GSH

☛ **CAS Number** :- 70-18-8

☛ **Detection Method** :- Colorimetric

☛ **Detection Range** :- 0.01 mmol/L – 0.5 mmol/L

☛ **Storage/Stability** :- Shipped and store at 4°C for 6 months.

☛ **Sample Types** :- Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

**Date Created**

2024/07/03