

ELISA Kit for Interleukin 8 (IL8)

Description

Catalog No.	SEA080Mu
Unit Size	48T, 96T, 96T*5, 96T*10, 96T*100
Product Category	ELISA (Quantitative)
Reactivity	Mus musculus (Mouse)
Sample	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and biological fluids
Assay Method	Double-antibody Sandwich
Assay Duration	3h
Sensitivity	The minimum detectable dose of this kit is typically less than 6.3pg/mL
Standard Curve Range	15.6-1,000pg/mL
Storage/Storage	-20°C/1 year; 4°C/6 months

Date Created

2024/07/03

LINK BIOTECH