

HRP-conjugated Affinipure Goat Anti-Mouse IgG (H+L)

Description

PRODUCT INFORMATION:-

☛ **Product Name** :- HRP-conjugated Affinipure Goat Anti-Mouse IgG (H+L)

☛ **Catalogue No.** :- FNSA-0003

☛ **Size** :- 100 μ l/500 μ l

☛ **Form** :- liquid

☛ **Conjugation** :- HRP

☛ **Stabilizer** :- 3mg/ml BSA (IgG-Free, Protease-Free)

☛ **Storage** :- 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.
-20 $^{\circ}$ C for 24 months , avoid repeated freezing and thawing.

☛ **Tested Application** :- ELISA, WB, IHC

☛ **Suggested Dilution** :- 1:1000-1:10000 for western blotting with ECL substrates
1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
1:20-1:50 for IHC

Date Created

2024/07/03