

## HRP-conjugated Goat Anti-Rabbit IgG

### Description

Alternative Names:	Â HRP-conjugated Goat Anti-Rabbit IgG
Catalogue No.	300SA
Size	1 ml
Applications	WB,IHC,ELISA WB 1:5000-1:10000 IHC 1:100-1:200
Recommended dilutions	ELISA Use at an assay dependent concentration NOTE: Optimal dilutions should be determined by the end user
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at 4°C. <a href="https://www.elkbiotech.com">https://www.elkbiotech.com</a> T:86-27-59760950 E: elkbio@elkbiotech.com
Note	This product is for research use only, not for diagnostic or therapeutic use.
<b>Date Created</b>	
2024/07/03	

LINK BIOTECH