

HRP-conjugated Rabbit Anti-Goat IgG

Description

Alternative Names: HRP-conjugated Rabbit Anti-Goat IgG

Catalogue No. 400SA

Size 1 ml

Applications WB,IHC,ELISA

WB 1:5000-1:10000

IHC 1:100-1:200

Recommended dilutions

ELISA Use at an assay dependent concentration

NOTE: Optimal dilutions should be determined by the end user.

Storage

PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. <https://www.elkbiotech.com> T:86-27-59760950 E: elkbio@elkbiotech.com

Note

This product is for research use only, not for diagnostic or therapeutic use.

Date Created

2024/07/03