

(R)-1-Isothiocyanato-4-(methylsulfinyl)butane

Description

Synonyms	L-Sulforaphane; (-)-Sulforaphane; (L)-Sulforaphane
Catalogue No.	BS631673
CAS	142825-10-3
Purity	98%
Molecular Formula	$C_6H_{11}NOS_2$
Molecular Weight	177.29
Storage	Inert atmosphere, Store in freezer, under -20°C
Pack Size	10mg 25mg 50mg 100mg 250mg
Documents	COA with the material

Date Created

2024/07/04

LINK BIOTECH