

Single Channel Variable Volume Micropipettes

Description

- Easily Modify Volume using Plunger
- The Plunger's superior spring mechanism ensures smooth, snag-free operation thanks to its thoughtful design. To easily change the instrument's loudness, just turn it.
- Utilize a Universal Tipcone with a Variety of Tips
- A Universal Tipcone improves the instrument's interoperability and makes it simple for it to operate with the majority of globally recognized standard tips.
- Perfectly adjust the volume.
- The ideal volume level is ensured and any unintentional adjustments are prevented with a gentle click sound at each increment. It also makes operating with one hand easier.
- Make Sure to Store with a Holder
- The instrument comes with a Holder that makes storage simple, effective, and secure.
- Reject Advice Using a Tip Ejector Made Easy
- Comfortable access and simple tip ejection are made possible with an integrated, streamlined tip ejector.

MODEL NO.	VOL. RANGE	INCREMENT (µL)	ACCURACY		CV	
			±%	±µL	±%	±µL
BSSLP-2	0.2-2.0 µL	0.01 µL	2	0.04	1.2	0.024
BSSLP-10	0.5-10 µL	0.02	1	0.1	0.5	0.05
BSSLP-20	2-20 µL	0.02	0.8	0.16	0.4	0.08
BSSLP-50	5-50 µL	0.1	0.8	0.4	0.4	0.2
BSSLP-100	10-100 µL	0.2	0.6	0.6	0.2	0.2
BSSLP-200	20-200 µL	0.2	0.6	1.2	0.2	0.4
BSSLP-1000	100-1000 µL	1.0	0.6	6	0.2	2
BSSLP-5000	0.5-5 ml	10.0	0.6	30	0.2	10
BSSLP-10000	1-10 ml	20.0	0.6	60	0.2	20

Date Created

2024/07/03