

Sterile Serological Pipettes

Description

Sufficient inventory, Short Delivery Time, New Price, High Cost Performance!

Fu Dau is one of the professional Serological Pipettes suppliers in China. Buy or wholesale high-quality products at competitive prices and get quotes from our factory.

Application scope: Accurately measure or transfer a fixed volume of liquid, widely used in cell culture, bacteriology, clinical medicine, scientific research and experiment fields.

Product material: USP-Class VI grade polystyrene raw material, produced in a GMP 100,000 grade clean workshop.

Product specifications: 1ml, 2ml, 5ml, 10ml, 25ml, 50ml, 100ml

Features:

- The tube is transparent and the scale is clear, which is easy to observe and read.
- Negative scale design provides a certain amount of additional sampling capacity.
- The end of the tube is marked with a color ring, which is convenient for quickly identifying the capacity.
- No pyrogen, no endotoxin, no cytotoxicity.
- The front and back double scale design is convenient for the operation of adding and subtracting samples.
- Irradiation sterilization, SAL 10⁻⁶.
- Smooth tube wall reduces liquid adhesion and improves sampling accuracy.
- Packing method: independent paper and plastic or in bulk.

Product Number	Volume	Sterilization	Packaging
SP010001	1 ml	Aseptic packaging	50 pieces
SP010002	2 ml	Aseptic packaging	50 pieces
SP010005	5 ml	Aseptic packaging	50 pieces
SP010010	10 ml	Aseptic packaging	50 pieces
SP010025	25 ml	Aseptic packaging	50 pieces
SP010050	50 ml	Aseptic packaging	50 pieces

Date Created
2024/07/03